

AMP USR *mini* - Portable Urine Strip Reader



Reliable

Urine Diagnostics

at Your Hand



AMP USB *mini*

Portable Urine Strip Reader

AMP USB mini is a modern urine strip reader, that enables instrumental reading of urine test strips literally everywhere. Dimensions and weight similar to a smart phone make the device really portable. The integrated high capacity Li-battery allows independent operation. Test results are saved in memory and can be printed with a wireless printer or downloaded via the USB port.



Technical Specifications

Measuring principle:	Light reflectance (LED)
Test parameters:	10 or 11
Throughput:	50 strips / hour normal mode 140 strips / hour fast mode
Display:	2.4" LCD screen
Operation:	5 touch keys
Printer:	wireless thermoprinter (optional)
Memory:	1.000 results
Dimensions, weight:	110 x 68 x 27 mm, 160 gr.
Power adapter:	100 - 240 VAC, 50/60 Hz
Power consumption:	less than 2 W
Battery capacity:	min. 6 hours operation



AMP Urine Test Strips and Quality Control Material

AMP Urine Strips US-10:

Leukocytes, Nitrite, Urobilinogen, Protein, pH, Blood, Specific Gravity, Ketones, Bilirubin, Glucose

AMP Urine Strips US-11:

same as US-10 plus Ascorbic Acid (Vitamine C)

AMP UroTrol:

Quality control material supplied in 10mL dropper bottles

AMEDA Labordiagnostik GmbH

Krenngasse 12

8010 Graz, Austria

Phone +43 316 69 80 69

Fax +43 316 69 80 69 12

E-mail office.graz@amp-med.com

www.amp-diagnostics.com



AMP Polska sp. z o.o.
AMP Diagnostics s.r.l.

Gdynska 31, 31-323 Kraków, Poland, Phone: +48 12 294 02 78, E-mail: office.krakow@amp-med.com
Str. Fabricii 118-A, 400632 Cluj-Napoca, Romania, Phone: +40 264 425 934, E-mail: office.cluj@amp-med.com